

focus on Chromatography

Analysis of Monoclonal Antibody Charge Variants by Capillary Zone Electrophoresis

Marcia R Santos, Global Tactical Marketing, Beckman Coulter Life Sciences, Brea, CA USA

Charge heterogeneity analysis is important in the characterisation of monoclonal antibodies because it provides important information about product quality and stability. Heterogeneity can be caused by such molecular adaptations as C-terminal lysine modification, deamidation, and post translational modification. One method for separating charge variants is capillary isoelectric focusing (cIEF) which provides information regarding isoelectric point variation for related molecular isoforms. cIEF requires electrophoretic separation in a coated capillary which helps suppress electro-osmotic flow (EOF) and prevents surface protein adsorption but can result in a fairly long cycle time and complex sample preparation.

We demonstrate here that by using a simple separation buffer system and a bare fused silica capillary, it is possible to obtain a highly resolved, reproducible separation of a representative monoclonal antibody in less than 12 minutes. We also demonstrate that while high resolution can be achieved using a short effective length (20cm), the fine structure of the monoclonal antibody used in this study can be revealed by increasing the effective length to 40cm.

Experimental

Capillary Zone Electrophoresis

Preparation of Separation Buffer: 0.05% HPMC, 380 mM EACA, 1.9 mM TETA

Preparation of 400mM EACA, 2mM TETA pH 5.7 Solution

5.25g of EACA (Sigma, cat. no. A-7824) and 30 μ l of TETA (triethylenetetramine) (Sigma, cat. no. 90460) were dispensed into a 100mL beaker containing 95mL of double distilled and deionised (ddi) water. The contents were allowed to dissolve completely. The pH was adjusted with glacial acetic acid to pH 5.7 + 0.05. The contents were quantitatively transferred to a 100mL volumetric flask and the volume was brought to 100mL with ddi water. The resulting solution was filtered through a 0.2 μ m filter.

Preparation of 1% HPMC Solution

One gram of HPMC (Sigma, cat. no. H-7509) was dispensed into a large beaker containing 100mL of ddi water. The solid was allowed to dissolve to completion overnight at room temperature.

Preparation of the Separation Buffer

Transfer 8.55mL of the 400mM EACA, 2mM TETA pH 5.7 buffer to a 15mL conical tube and add 450 μ L of 1% HPMC solution. Mix well. Prepare fresh before each use. The amount prepared is enough for 9 separations.

Capillary Wash Solution

0.1 N HCl solution (Fluka, part no. 94015).

Sample

MAB X is a representative, therapeutic grade monoclonal antibody. A 20mg/mL stock solution of MAB X was diluted to a final concentration of 1mg/mL in ddi water and used in both cIEF and CZE separations.

Capillary

Type: Bare Fused Silica; 50 μ m i.d. x 360 μ m o.d.

Effective length: 40cm, total length 50cm for high resolution method and 20cm effective length, total length 30cm for the fast separation method.

Instrument

PA 800 plus Pharmaceutical Analysis System equipped with either a UV or PDA detector was used in these experiments. UV detection was configured with a 214nm filter, with a data rate set at 4Hz to perform the CZE experiment and 280nm filter to perform cIEF. PDA detection was also used in duplicate CZE experiments for later comparison with the UV/vis data. The sample storage and cartridge temperature were kept at 20 and 25°C respectively.

The capillary conditioning method was performed using a series of rinses followed by sample injection and separation as follows: 50 psi rinse of 0.1 N HCl solution for 5 min followed by buffer rinse at 50 psi for 5 min. The sample was injected hydrodynamically at 0.5 psi for 10 seconds. The separation was performed using normal polarity at 30 kV for 30 min.

There are two distinct separation methods for high speed and high resolution analysis. The high speed method used for a 30 cm bare fused silica capillary consists of 50 psi rinse of 0.1 N HCl solution for 5 min followed by buffer rinse at 50 psi for 5 min. The sample is then injected hydrodynamically at 0.5 psi for 10 seconds. The separation was performed using normal polarity at 30 kV. For the molecule used in this study the duration of the separation was only 12 min. The high resolution method used with a 50cm bare fused silica capillary is very similar to the high speed except the buffer rinse step is 10 min long and the separation step is 40 min long. The separation time may vary depending on the pI (protein isoelectric point) of the MAB under test.

The shutdown method consists of a 50 psi rinse with 0.1 N HCl solution for 5 minutes, followed by a wait step where the ends of the capillary are immersed in ddi water and finally a 'lamp off' step.

Results And Discussion

Charge heterogeneity analysis of MAB X was performed using both cIEF and CZE. A typical cIEF separation of MAB X illustrates a complex profile with charge isoforms ranging from pI 6.3 to pI 6.91 (Figure 1). It wasn't clear whether this MAB would separate efficiently using CZE since separation buffer was pH 5.7, potentially leading to solubility issues as discussed by He et al[1].




Figure 1. cIEF profile of MAB X. MAB concentration was 1mg/mL.

Even though the protein is somewhat acidic and with less than one pH unit difference from the pH of the separation buffer, it was possible to obtain a complete high resolution separation in only 12 minutes (Figure 2).




Figure 2. CZE profile and peak integration strategy for MAb X. Peaks were grouped as either Basic, Main or Acidic variants. CZE conditions: MAb X 1mg/mL; Separation buffer: 0.05% HPMC, 380 mM EACA (epsilon amino caproic acid), 1.9 mM TETA; bare fused silica capillary, 30 kV separation voltage.

Results illustrate that high resolution separation can be achieved using the CZE methodology (Figure 3). In addition, the resolution in the basic region surpasses that of cIEF. An important attribute of CZE is not only the high resolution of this separation but also the fast separation time. A complete charge heterogeneity profile can be obtained in 12 minutes using CZE compared to 23 minutes required by cIEF.




Figure 3. Comparison between cIEF and CZE of MAb X.

Given the importance of assay repeatability, a number of CZE separations were performed sequentially to illustrate reproducibility (Figure 4). Data integration on these data was performed so that basic, main, and acidic peak groups could be designated as shown in Figure 2.




Figure 4. CZE separation of MAb X. Nine consecutive separations of MAb X were performed using UV detection. MAb X concentration was 1mg/mL; separation buffer: 0.05% HPMC, 380 mM EACA, 1.9mM TETA; bare fused silica capillary, 30 kV separation voltage.

This CZE separation resulted in exceptional reproducibility for a number of consecutive runs; better than 3% RSD for corrected peak areas and better than 0.55% RSD for migration time for each of the basic, main and acidic peak groups (Table 1). Corrected area is commonly used in capillary electrophoresis because each peak travels at different velocity past the detection point. Here the corrected area was calculated by dividing the area under the peak by its velocity.

Table 1. Raw data for CZE separation of MAb X. Reproducibility of the CZE separation for corrected area and migration time for basic, acidic and main grouped peaks.

CZE of Mab X	Corrected Area			Migration Time		
	Basic	Main	Acidic	Basic	Main	Acidic
Mab X run #1	12899.79	21950.23	4363.77	8.24	8.99	9.81
Mab X run #2	12382.87	21373.5	4155.31	8.13	8.87	9.67
Mab X run #3	12572.88	21441.31	4333.81	8.18	8.91	9.73
Mab X run #4	12483.24	21438.48	4383.71	8.22	8.97	9.78
Mab X run #5	13410.81	22744.52	4581.26	8.15	8.89	9.7
Mab X run #6	12900.86	21703.76	4418.44	8.23	8.98	9.8
Mab X run #7	12906.9	21945.95	4437.96	8.23	8.98	9.8
Mab X run #8	13118.91	22094.07	4569.69	8.23	8.98	9.8
Mab X run #9	13046.59	21908.14	4501.46	8.24	8.99	9.82
Mean:	12858.09	21844.44	4416.16	8.21	8.95	9.77
Std Dev:	328.39	428.21	130.96	0.04	0.05	0.05
%RSD:	2.55	1.96	2.97	0.51	0.53	0.55

The CZE method also has the flexibility to offer the same performance regardless of the type of detection. Comparison between two CZE separations of MAb X using a photodiode array detector and UV detector illustrates the same resolution can be obtained with both detectors. Additionally, the area % composition is preserved across both detectors (Figure 5).




Figure 5. Comparison of UV and PDA in CZE separation of MAb X. CZE was performed using either UV detection (red trace) or PDA detection (black trace). Separation conditions are the same as those described in Figure 1.




Figure 6. Comparison between 20 and 40cm capillary effective length. CZE separations were performed as described previously using either 40cm capillary (bottom trace) or a 20cm capillary (inset).

In cases where resolution is more important than speed of separation, by increasing the effective length of the capillary, it is possible to unveil fine structures that were not obvious in the fast separation method. As shown in *Figure 6*, the separation of MAb X using a capillary with 40cm effective length, the pay off of a high resolution is clear even though the total separation time increased significantly.

The increase in resolution may be a benefit when analysing a degraded form of MAb X, exposed to 60°C for a period of 48 hours. However, only by increasing the effective length was it possible to see the fine structure of acidic peaks and basic peaks formed due to temperature instability (*Figure 7*).




Figure 7. MAb X degraded by exposure to high temperature (60°C) for 5 days. CZE separations were performed as described previously using either 40cm capillary (bottom trace) or a 20cm capillary (inset).

Conclusion

Charge heterogeneity of protein isoforms can easily be determined using CZE and cIEF. While cIEF is commonly used to obtain charge heterogeneity and pI information for heterogeneous isoforms, we illustrate here that CZE is capable of providing faster, highly reproducible separations for a representative monoclonal antibody.

This CZE method can be optimised by increasing the effective length of the capillary to increase peak resolution. Quality of this separation was not affected by detection method and reproducibility of this CZE assay is excellent across basic, main, and acidic species. Given the increasing need for sample throughput, separation quality, and repeatability, CZE can be positioned as a faster alternative to cIEF.

References

1. Yan He, Colleen Isele, Weiyang Hou, Margaret Ruesch. 2011. Rapid analysis of charge variants of monoclonal antibodies with capillary zone electrophoresis in dynamically coated fused-silica capillary. *Journal of Separation Science*, vol. 34: pages 548-555.
2. Ingrid D. Cruzado-Park, Scott Mack and Chitra K. Ratnayake, A Robust cIEF Method: Intermediate Precision for the pH 5-7 Range, PN A-12015

Acknowledgment

The author would like to thank Yan He from Pfizer and Bernd Moritz from Hoffman LaRoche for valuable discussions.

Read, Print, Share or Comment on this Article at: Labmate-Online.com/Articles



Online Filtration Capability Makes HT-GPC Safer

According to Dr Mark Jordi, President at Jordi Labs LLC, the online filtration capability offered by Viscotek high temperature gel permeation chromatography (HT-GPC) systems from **Malvern Instruments** has removed the need for dangerous and difficult manual filtration during high temperature polymer analysis.

"GPC is an integral part of Jordi Labs' analytical capabilities," said Dr Jordi. "The Viscotek systems provide a richness of information not obtained with other instruments. The online filter for the high temperature GPC instruments works well. We previously performed this manually and it was both dangerous and difficult. Online high temperature filtration is a big improvement."

Dr Jordi continued; "As experts in polymer characterisation, offering over 30 years of experience in the field, Jordi Labs use Viscotek GPC systems for routine quality control, failure analysis, and product deformation. The ability to compare polymer molecular weight, intrinsic viscosity and molecular size in a single measurement is powerful for identifying subtle differences in a polymer system. We especially like the ability to prepare and evaluate the Mark Houwink plots to compare branching in polymer systems. Two high temperature systems are also used extensively for polyethylene and polypropylene characterisation."

Jordi Labs specialises in performing laboratory analysis using High Temperature GPC, also known as High Temperature Gel Permeation Chromatography.

[MORE INFO.](#) 31



Automation-Friendly Tabless Tube Holder

Phenomenex Inc introduces two new Tabless Tube Holder accessories for its Strata® and Strata™-X lines of silica-based and polymeric solid phase extraction (SPE) sorbents. These accessories hold SPE tubes in 96-well spacing to integrate with automated systems.

All Strata and Strata-X products are now available in the tabless 1mL tube format, enabling the user to arrange multiple sorbents within the holders, adding flexibility for method development. Users can also remove and replace a single SPE tube should they make an error, which cannot be done with standard 96-well SPE plates. The Tabless Tube Holders are compatible with the Phenomenex vacuum manifold as well as positive pressure systems.

"Our new Tabless Tube Holders are easy to use and add flexibility to the workflow options using our Strata and Strata-X sorbents," commented Erica Pike, Brand Manager for Phenomenex. "Our growing line of Strata and Strata-X products meets the sample preparation needs of a wide range of applications, including complex biological matrices."

Strata and Strata-X SPE sorbents simplify the method development process for fast and efficient sample preparation prior to chromatography. They remove unwanted contaminants, including phospholipids, and are offered in a number of unique selectivities to cover a diverse range of analytes.



[MORE INFO.](#) 32

HPLC Assay of Gentamicin Sulphate Congeners in Just 15 Minutes

Thermo Fisher Scientific is pleased to announce a direct-detection reversed-phase high-performance liquid chromatography chromatographic assay method that separates the congeners in gentamicin sulphate and gentamicin sulphate-containing ointments, solutions, and creams in just 15 minutes.

Application Note 1005: Gentamicin Sulfate Assay by HPLC with Charged Aerosol Detection demonstrates that this approach is faster than previously published methods, with excellent separation between the gentamicin congeners and resolution of many impurities that can be present in fermentation-manufactured products.

This method uses a silica-based polar-embedded column that provides enhanced hydrolytic stability from pH 1.5–10 and allows resolution of a wide variety of polar and nonpolar analytes. Use of this column allows a reduced concentration of trifluoroacetic acid, thus improving method ruggedness.

Use of charged aerosol detection eliminates the need for derivatisation. With three sample matrices analysed, the quantification of gentamicin base had good recoveries and was consistent with the label claim of these products, both indicators of an accurate assay.

[MORE INFO.](#) 33